

How to Calculate a pH Probe Slope Percentage

Calculating a pH slope percentage will allow you to verify that your pH probe is functioning properly.

Often, you will obtain a straight calibration curve on a pH probe and obtain consistent readings, but these readings are erroneous.

This is because the Offset of the curve is outside of the range of acceptable values.

Generally a slope between 90 and 105% is acceptable.

Outside of this range, you may need to replace your pH probe.

To calculate the Slope Percentage of your pH Probe:

- 1) Read the mV potential generated by the electrode in two separate calibration buffers.
 - 2) Calculate the difference in mV's between the higher buffer value and the lower buffer value.
 - 3) Divide the difference in mV potential by the difference in pH unit
(**EXAMPLE:** the difference in pH unit between 7.01 buffer and 4.01 buffer is $7.01 - 4.01 = 3$).
 - 4) Divide this number by the theoretical maximum (59.16 mV/pH unit @ 25 °C) and multiply by 100 to get a percentage value.
-

Sample Calculation:

pH electrode generated -15 mV in pH 7.01 buffer and +160 mV in pH 4.01 buffer.

$$160 \text{ mV} - (-15 \text{ mV}) = 175 \text{ mV}$$

$$175 \text{ mV} / 3 = 58.33 \text{ mV/pH unit}$$

$$(58.33 / 59.16) \times 100 = 98.6\% \text{ slope. } \textit{\textbf{This is your pH slope percentage.}}$$