

**Above and Beyond with Carolina's
AP® Biology Series: Explore the
Options!**

Andrew Uy
Carolina Teaching Partner

**NSTA 2009
New Orleans, LA**



Questions to Focus Session

- **What is your level of knowledge for AP[®] Biology?**
- **Do you meet the needs of your students through laboratory investigations?**
- **What topics interest you most in AP[®] Biology?**

Objectives

- **To introduce Carolina's newly revised and improved line of AP[®] Biology lab kits**
- **To demonstrate the wide variety of AP[®] Biology topics offered for laboratory investigation**
- **To use materials and procedures from 5 different Carolina[™] AP[®] Biology kits**

Carolina's AP[®] Environmental Kits

12 Labs Revised and Improved!

- **Diffusion and Osmosis**
- **Enzyme Catalysis**
- **Mitosis and Meiosis**
- **Plant Pigments and Photosynthesis**
- **Cell Respiration**
- **Genetics of *Drosophila***
- **Population Genetics and Evolution**
- **Transpiration**
- **Physiology of the Circulatory System**
- **Animal Behavior**
- **Dissolved Oxygen and Aquatic Primary Productivity**

Jigsaw Activities

Group A:

Observing Catalase Reaction

Group B:

Leaf Pigment Chromatography

Group C:

Cell Respiration

Group D:

Sexing *Drosophila* and Examining Phenotypes

Group E:

Heart Rate of *Daphnia*

Safety Issues with Materials

- **Most procedures suggest the use of safety glasses, aprons, and gloves**
 - **Acids, stains in use**
- **Disposal**
 - **Follow instructions in handout for proper disposal procedure for your activity**



Reporting Out: Attendee Ideas

- **AP® Biology topics covered**
- **How materials met needs**
- **How would use in own classroom**
- **Value to students**
- **Need for classroom**


Group A: Observing Catalase

Activity Objective


To observe the action of an enzyme and the effects of heat on enzyme activity

Adapted from Lab 2:
Enzyme Catalysis

1. Observe the reactions of catalase on H_2O_2



2. Does living tissue contain H_2O_2 ?



Group B: Leaf Pigment Chromatography


Activity Objective

To separate plant pigments using paper chromatography and calculate R_f values from data.


Adapted from Lab 4:
*Plant Pigments and
Photosynthesis*

1. Prepare pigment origin

2. Place paper in solvent jar




pencil marks
1.5 cm from
edge



Group B: Leaf Pigment Chromatography

3. Measure migration distances and calculate R_f values



Sample data table

	Distance from origin (mm)		
Solvent Front	68		
Band #		Band Color/Identification	R_f Value
1	8	<i>dark green</i>	.12
2	13	<i>green</i>	.19
3	20	<i>yellow</i>	.29
4	65	<i>yellow-orange</i>	.96

$$R_f = \frac{\text{distance of pigment from origin}}{\text{distance of solvent front from origin}}$$


Group C: Cell Respiration

Activity Objective

To set up a respirometer to measure the consumption of oxygen by germinating peas

1. Respirometer setup

germinating peas →
nonabsorbent cotton →
absorbent cotton w/ KOH →



Adapted from Lab 5:
Cell Respiration

← graduated pipet

Group C: Cell Respiration

2. Place respirometer in water bath



glass beads only

dry peas + glass beads

germinating peas

Group C: Cell Respiration

- As oxygen is consumed by seeds, water enters the pipet. Measure decrease in oxygen volume over time.



Sample Table 1: Respiration of Peas at Room Temperature


°C	Time (Min)	Respirometer 1 Germinating Peas			Respirometer 2 Dry Peas + Beads			Respirometer 3 Beads Only	
		V of Pipet	ΔV	Corrected ΔV	V of Pipet	ΔV	Corrected ΔV	V of Pipet	ΔV
24	0	0.69	–	–	0.72	–	–	0.74	–
24	5	0.60	0.09	0.06	0.71	0.02	0.00	0.72	0.02
24	10	0.41	0.28	0.24	0.69	0.07	0.03	0.70	0.04
24	15	0.30	0.39	0.34	0.67	0.10	0.05	0.69	0.05
23	20	0.18	0.51	0.46	0.65	0.16	0.11	0.69	0.05

$\Delta V = V$ at Time 0 – V at time of current reading

Corrected $\Delta V = \Delta V$ (for Respirometer 1 or Respirometer 2) – ΔV of Respirometer 3

Group C: Cell Respiration

4. Graph of sample data




Group D: Examining *Drosophila*

Activity Objective

To determine the sex of fruit flies and observe contrasting phenotypes

Adapted from Lab 5:
Genetics of Drosophila

1. Distinguish males from females



Group D: Examining *Drosophila*

2. Observe contrasting phenotypes



wild type



A. sepia



B. white



C. vestigial

Phenotype	Phenotype
Wild Type: red eyes	Mutant A: brown eyes
Wild Type: red eyes	Mutant B: white eyes
Wild Type: normal wings	Mutant C: short wings




Group E: Heart Rate of *Daphnia*

Activity Objective

To use *Daphnia* to study the effect of environmental temperature on heart rate

Adapted from Lab 10:
*Physiology of the
Circulatory System*

Double concave slide assembly



Group E: Heart Rate of *Daphnia*


Sample Data Table

Temperature °C	Heartbeats/10 sec	Heart rate in beats/min (Heartbeats/10 sec × 6)
7	25	150
12	38	228
15	53	318
18	61	366
24	72	432

Group E: Heart Rate of *Daphnia*

Sample Graph

Daphnia Heart Rate at Different Temperatures





What Kits Have You Seen Today?

- Enzyme Catalysis
- Plant Pigments and Photosynthesis
- Cell Respiration
- Genetics of *Drosophila*
- Physiology of the Circulatory System



Other Quality AP[®] Bio Kits Available



- Diffusion and Osmosis
- Mitosis and Meiosis
- Molecular Biology
- Population Genetics and Evolution
- Transpiration
- Animal Behavior
- Dissolved Oxygen and Primary Productivity





Summary

- **Carolina's kits offer a wide variety of AP[®] Biology topics**
- **Newly revised kits provide inquiry-based learning method**
- **Materials are great quality and easy to use**
- **?????**

Evaluations: Share Your Thoughts!

- **Scale = 1 to 10**
- **10 = Outstanding**
- **9 = Above Average**
- **8, 7 = Average**
- **6, 5, 4 = Below Average**
- **3, 2, 1 = Well Below Average**
- **Please provide comments!**





Carolina Biological Supply Company

**Thank you for investing your time in
our training program.**

**For all your classroom needs, check out
our website, www.carolina.com.**

Enjoy the rest of the conference!