





BAL	QTY	UM	P/N	DESCRIPTION
1	1	EA	2065900	CABINET 2946K 2946K
2	1	EA	2065939	APRON ASSY OAK 4 STUDENT LABVIEW
3	1	EA	101013	TOP PANEL 10050X868 BLK FTR LABVIEW
4	2	EA	138-W	LEG WING OAK 38
5	1	EA	2066100	ELEC ASSY OPT 13A 125V/BLK
6	2	EA	100170	BAG MINI GRIP 4X4 CLEAR
7	2	EA	5200092	TABLE LEG COMPONENTS
8	2	EA	2066620	LEG HARDWARE ASSY C/INTN BUTS NUTS
9	2	EA	2068210M-LH	DOOR HALL KNC BR 44 52X29 25 LH
10	1	EA	2068210M-RH	DOOR HALL KNC BR 44 52X29 25 RH
11	1	EA	2065912	HINGE ASSY FILLOVER/VA/ST
12	4	EA	2065907	CATCH ROLLER ASSY
13	2	EA	2065922	PULL ASSY DOOR 28MM W/HR CHROME
14	2	EA	100283	RUBBER BASE INCLUDING (SOLD BY FOOD CORNER BASE ASSY SST 400NERSV/SCW)
15	1	EA	100074	EXTDRE LEG -WSA -PW(MULTI-SERVICE)
16	2	EA	206581	SINK ASSY FPOXY STRIP STRIP TRAP
17	1	EA	100170	BAG MINI GRIP 4X4 CLEAR
18	2	EA	200397-2218	SHIELD ADJ K03/EHT 75X11.00X22.18
19	2	EA	100072	SHIELD CLAMP ASSY 30mm PIN
20	1	EA	219218	BRACKET HINGE PARSNT TO CABINET ASSY



ITEMS NOT VISIBLE (7, 20)



3D VIEW  
SCALE: 1" = 12"  
(shown for clarity)

NOTE:  
ITEMS 8-10 ARE USED TO ATTACH LEGS TO CABINET

REV	DESCRIPTION	DATE	BY	CHKD	APP'D
1	ADD ITEMS 2, 3, 20 FROM CABINET PART 23 - 08.18.18.10021	8/1/2014	JABR/IA		
2	REVISED DRAW SHEET WIDTH & LOCK PART NUMBER	6/27/04	STUSCHT		
3	ITEM 10 WAS QTY 1, ITEM 26 WAS 4, & ADDED 28	7-28-04	STUSCHT		
4	Removed back-roll and added corners	2-18-05	STUSCHT		

DIVERSIFIED WOODCRAFTS, INC  
LABVIEW 4-STUDENT WORKSTATION W/FIX  
2946K  
REV 4