




BAL	QTY	UM	P/N	DESCRIPTION
1	1	EA	181-4822	SINK BASE 48Wx22D
2	1	EA	110-3022	4 DRAWER/RH DOOR 30Wx22D
3	1	EA	207824	TOP, EPOXY, 1.00X108X24, BLK, FIXT, 4PC
4	1	EA	100699	SINK, EPOXY, 130, BLK
5	1	EA	204561	SINK, ASSY, EPOXY, STRNR, STPPR, TRAP
6	2	EA	100074	FIXTURE, L6S-W5A-DIV (MULTI-SERVICE)
7	12,292	FT	100283	RUBBER BASE MOLDING (SOLD BY FOOT)
8	1	EA	204513	CORNER BASE ASSY, SST, 2CORNERSW/SCRW
9	1	EA	116-3022	4 DRAWER/LH DOOR 30Wx22D
10	2	EA	100033	ELECTRICAL BOX 3.00X2X2.5 DEEP
11	2	EA	100034	ELEC. OUTLET, GR, 20A 125V, BLK



NOT SHOWN

⑩ ⑪



REV	DESCRIPTION	DATE	APPROVED	DESIGNER	CHECKER	DATE
				STUSCHY		10/12/2004

DIVERSIFIED WOODCRAFTS, INC	
WALL SERVICE BENCH 108Wx24DX34H W/CABINET 181-4822, 110-3022 & 116-3022	
DRAWING NUMBER	REVISION
3236K	
SCALE 1"=8"	SHEET 1 OF 1

THIS DRAWING IS THE SOLE PROPERTY OF DIVERSIFIED WOODCRAFTS. IF ANY NOT BE ALTERED OR CHANGED WITHOUT PERMISSION