

Prepared Slides

Care Instructions

Proper storage and handling are essential for preserving the quality of microscope slides. Please follow the simple guidelines below when storing your slides:

- Store slides in a container made for holding microscope slides
- Keep slides away from heat or bright light
- Store slides in the horizontal (flat) position. If stored vertically (on edge), the coverslip or specimen may drift
- If you have a slide set, store the set box in the upright (vertical) position, which will allow your slides to lay flat when not in use.
- Do not stack slides on top of one another or apply pressure to the coverslips

Problems? We hope not, but if so contact us. We want you to have a good experience.

Orders and replacements: 1-800-334-5551 then select Customer Service

Technical Support and Questions: customer_service@carolina.com



© 2009 Carolina Biological Supply Company